

(R)-(+)-Atenolol


Catalog No: #S0007L

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	(R)-(+)-Atenolol
Brief Description	Inhibitors
Purification	99.88%
Target Name	Adrenergic Receptor agonist
Calculated MW	266.34
Formulation	C14H22N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: CC(C)NC[C@H](COC1=CC=C(C=C1)CC(=O)N)OPathways: GPCR/G Protein Receptor: α -adrenergic receptor Boiling pt: 508C Melting pt: 148-152C Solubility: soluble in Water Appearance: White

Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Nuttall SL, et al. J Clin Pharm Ther. 2003 Jun;28(3):179-86.

Note: This product is for in vitro research use only