

Phenethyl alcohol

Catalog No: #S0011


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Phenethyl alcohol |
| Brief Description | Inhibitors |
| Purification | 99.82% |
| Target Name | Others |
| Calculated MW | 122.16 |
| Formulation | C8H10O |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:ImmunologySMILES:OCCC1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:219-221CMelting pt:-27CSolubility:soluble in WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lingappa BT, et al. Science. 1969 Jan 10;163(3863):192-4.

Note: This product is for in vitro research use only