

Phenethyl alcohol


Catalog No: #S0011

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phenethyl alcohol
Brief Description	Inhibitors
Purification	99.82%
Target Name	Others
Calculated MW	122.16
Formulation	C8H10O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OCCC1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 219-221°C
Melting pt: -27°C
Solubility: soluble in Water
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lingappa BT, et al. Science. 1969 Jan 10;163(3863):192-4.

Note: This product is for in vitro research use only