

Mecarbinat

Catalog No: #S0013


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mecarbinat
Brief Description	Inhibitors
Purification	99.12%
Target Name	Others
Calculated MW	233.27
Formulation	C ₁₃ H ₁₅ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCOC(=O)C1=C(C)N(C)C2=C1C=C(O)C=C2
 Pathways: Others
 Receptor: Others
 Boiling pt: 399.8°C
 Melting pt: 208-212°C
 Solubility: DMSO: 47 mg/mL (201.49 mM)
 Appearance: Gray Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Shadurski? KS, et al. Farmakol Toksikol. 1983 Mar-Apr;46(2):115-20.

Note: This product is for in vitro research use only