

Mephenesin


Catalog No: #S0015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Mephenesin |
| Brief Description | Inhibitors |
| Purification | 99.81% |
| Target Name | Others antagonist |
| Calculated MW | 182.22 |
| Formulation | C10H14O3 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:OtherSMILES:CC1=CC=CC=C1OCC(O)COPathways:OthersReceptor:OthersBoiling pt:153-154CMelting pt:69-71CSolubility:DMSO: 10 mMAppearance:White Cotton PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Prescire Int. 2013 May;22(138):127-8.

Note: This product is for in vitro research use only