

## Mexenone


Catalog No: #S0016

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mexenone
Brief Description	Inhibitors
Purification	99.40%
Target Name	Others
Calculated MW	242.28
Formulation	C15H14O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: COC1=CC(O)=C(C=C1)C(=O)C1=CC=C(C)C=C1 Pathways: Others Receptor: Others Boiling pt: 400.1°C Melting pt: 99-102°C Solubility: Water: 10 mM Appearance: Faint Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Torres V, Correia T. Contact Dermatitis. 1991 Aug;25(2):126-7.

Note: This product is for in vitro research use only