

## Phenolphthalein

Catalog No: #S0017


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Phenolphthalein
Brief Description	Inhibitors
Purification	99.09%
Target Name	Others
Calculated MW	318.33
Formulation	C20H14O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: OC1=CC=C(C=C1)C1(OC(=O)C2=CC=CC=C2)C1=CC=C(O)C=C1 Pathways: Others Receptor: Others Boiling pt: 548.7°C Melting pt: 258-263°C Solubility: Water: 400 mg/L Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Uesawa Y, et al. Arch Toxicol. 2007 Mar;81(3):163-8.

Note: This product is for in vitro research use only