

## Vanatone

Catalog No: #S0021


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Vanatone
Brief Description	Inhibitors
Purification	99.38%
Target Name	Others
Calculated MW	150.17
Formulation	C9H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CCC(=O)C1=CC=C(O)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: 152-154°C  
Melting pt: 36-38°C  
Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Cizmrikov R, Kozlovsk? J. Ceska Slov Farm. 1994 Oct;43(5):226-8.

Note: This product is for in vitro research use only