

Vanatone

Catalog No: #S0021


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vanatone
Brief Description	Inhibitors
Purification	99.38%
Target Name	Others
Calculated MW	150.17
Formulation	C9H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCC(=O)C1=CC=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: 152-154°C
Melting pt: 36-38°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Cizmrikov R, Kozlovsk? J. Ceska Slov Farm. 1994 Oct;43(5):226-8.

Note: This product is for in vitro research use only