

Taurine


Catalog No: #S0022

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Taurine |
| Brief Description | Inhibitors |
| Purification | 99.92% |
| Target Name | LDL inhibitor |
| Calculated MW | 125.15 |
| Formulation | C2H7NO3S |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Immunology
SMILES: NCCS(=O)(=O)O
Pathways: Metabolism
Receptor: LDL; VLDL
Boiling pt: Melting pt: >300°C
Solubility: Water: 20 mg/mL (159.8 mM)
Appearance: White Cotton Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Huxtable R.J. Prog Neurobiol, 1989, 32(6), 471-533.

Note: This product is for in vitro research use only