

Sulfalen

Catalog No: #S0026


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfalen
Brief Description	Inhibitors
Purification	99.29%
Target Name	Autophagy inhibitor
Calculated MW	280.3
Formulation	C11H12N4O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:COC1=NC=CN=C1NS(=O)(=O)C1=CC=C(N)C=C1Pathways:AutophagyReceptor:DHPSBoiling pt:488.6CMelting pt:270CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brumfitt W, Hamilton-Miller JM. J ChemOthers. 1993 Dec;5(6):465-9

Note: This product is for in vitro research use only