

4-Amino-3-hydroxybutyric Acid


Catalog No: #S0028

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Amino-3-hydroxybutyric Acid
Brief Description	Inhibitors
Purification	99.87%
Target Name	GABA Receptor agonist; Others
Calculated MW	119.12
Formulation	C4H9NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: NCC(O)CC(=O)O
Pathways: Membrane transporter/Ion channel; Others; Neuroscience
Receptor: GABAR; GHB receptor
Boiling pt: 374.5°C
Melting pt: 223°C
Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Candela M, et al. J Chromatogr A. 2000 Aug 25;890(2):273-80.

Note: This product is for in vitro research use only