

Miltefosine

Catalog No: #S0033


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Miltefosine
Brief Description	Inhibitors
Purification	99.45%
Target Name	PKC inhibitor
Calculated MW	407.58
Formulation	C21H46NO4P
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CCCCCCCCCCCCCCCCOP(O)(=O)OCC[N+](C)(C)CPathways:Cytoskeletal SignalingReceptor:PKCBoling pt:486C?Melting pt:149-153CSolubility:Ethanol: 82 mg/mL (201.19 mM); Water: 82 mg/mL (201.19 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Uberall F, et al. Y Res, 1991, 51(3), 807-812.

Note: This product is for in vitro research use only