

## Cortisone acetate

Catalog No: #S0034


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cortisone acetate
Brief Description	Inhibitors
Purification	99.74%
Target Name	GR inhibitor
Calculated MW	402.49
Formulation	C23H30O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CC(=O)OCC(=O)[C@]1(CC[C@@H]2[C@@]1(CC(=O)[C@H]3[C@H]2CCCC4=CC(=O)CC[C@]34C)O)Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:237-240CSolubility:DMSO: 7 mg/mL (17.39 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Grossman R, et al. Ann N Y Acad Sci. 2006, 1071:410-21.

Note: This product is for in vitro research use only