

Dienestrol


Catalog No: #S0035

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dienestrol
Brief Description	Inhibitors
Purification	99.05%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	266.34
Formulation	C18H18O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: C/C=C(/C(=C/C)/C1=CC=C(C=C1)O)\C2=CC=C(C=C2)O
 Pathways: Endocrinology/Hormones
 Receptor: ER
 Boiling pt: 395.4°C
 Melting pt: 229°C
 Solubility: DMSO: 10 mM
 Appearance: Faint Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Juriansz RL, et al. Y Res. 1988 Jan 1;48(1):14-8.

Note: This product is for in vitro research use only