

Halazone


Catalog No: #S0037

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Halazone
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	270.09
Formulation	C7H5Cl2NO4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC(=O)C1=CC=C(C=C1)S(=O)(=O)N(Cl)Cl Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. HACKENBERG HW. Z Gesamte Hyg. 1964 Apr;10:241-51.

Note: This product is for in vitro research use only