

Quinestrol

Catalog No: #S0040


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Quinestrol
Brief Description	Inhibitors
Purification	99.55%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	364.53
Formulation	C25H32O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@]2(C#C)O)CCC4=C3C=CC(=C4)OC5CCCC5Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:502.5C?Melting pt:107 - 108CSolubility:Water: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shyu C, et al. Toxicol Appl Pharmacol. 2010 Nov 12.

Note: This product is for in vitro research use only