

Amidrine


Catalog No: #S0043

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amidrine
Brief Description	Inhibitors
Purification	99.10%
Target Name	Others
Calculated MW	129.24
Formulation	C8H19N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCCC(C)CCN
Pathways: Others
Receptor: Others
Boiling pt: 155°C
Melting pt: -100°C
Solubility: DMSO: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. SEMENZA F. Boll Soc Ital Biol Sper. 1951 Mar;27(3 bis):354-7.

Note: This product is for in vitro research use only