

## Diphenylpyraline hydrochloride

Catalog No: #S0059


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diphenylpyraline hydrochloride
Brief Description	Inhibitors
Purification	99.26%
Target Name	Histamine Receptor antagonist
Calculated MW	317.85
Formulation	C19H23NOHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: Cl.CN1CCC(CC1)OC(C1=CC=CC=C1)C1=CC=CC=C1  
Pathways: GPCR/G Protein; Immunology/Inflammation; Neuroscience  
Receptor: H1 receptor  
Boiling pt: Melting pt: 203-205°C  
Solubility: Soluble in DMSO  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Lapa GB, et al. Eur J Pharmacol. 2005, 506(3):237-40.

Note: This product is for in vitro research use only