

Adenine


Catalog No: #S0064

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenine
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	135.13
Formulation	C ₅ H ₅ N ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: NC1=NC=NC2=C1N=CN2
Pathways: Neuroscience
Receptor: Adenine phosphoribosyltransferase
Boiling pt: >360°C
Solubility: Ethanol:

References

1. Barrett C, et al. Assay Drug Dev Technol. 2006 Dec;4(6):661-9.

Note: This product is for in vitro research use only