

Dichlorophen

Catalog No: #S0066


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dichlorophen
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	269.13
Formulation	C13H10Cl2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC1=C(C(C2=C(O)C=CC(Cl)=C2)C=C(Cl)C=C1
Pathways: Others
Receptor: Others
Boiling pt: 168-172°C
Melting pt: 168-172°C
Solubility: Water: 30 mg/L
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Langrand J, et al. Clin Toxicol (Phila). 2013 Mar;51(3):178-81.

Note: This product is for in vitro research use only