

## Oxiniacetic Acid

Catalog No: #S0067


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Oxiniacetic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	139.11
Formulation	C6H5NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: OC(=O)C1=C[N+](O-)=CC=C1  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 469°C  
 Melting pt: 254-255°C  
 Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Atac A, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2012 Jan;85(1):145-54.

Note: This product is for in vitro research use only