

Uracil

Catalog No: #S0069


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Uracil
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	112.09
Formulation	C ₄ H ₄ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1NC=CC(=O)N1Pathways:OthersReceptor:RNABoling pt:Melting pt:>300CSolubility:Ethanol: 22 mg/mL (196.27 mM); DMSO: 22 mg/mL (196.27 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Imming P, et al. Nat Rev Drug Discov. 2006 Oct; 5(10):821-34.

Note: This product is for in vitro research use only