

## Homosalate

Catalog No: #S0071


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Homosalate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	262.35
Formulation	C16H22O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: CC1CC(CC(C)(C)C1)OC(=O)C1=CC=CC=C1O Pathways: Others Receptor: Others Boiling pt: 341.1°C Melting pt: Solubility: DMSO: 55mg/mL Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kim TH, et al. J Toxicol Environ Health A. 2014;77(4):202-13.

Note: This product is for in vitro research use only