

Galanthamine HBr

Catalog No: #S0086


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galanthamine HBr
Brief Description	Inhibitors
Purification	99.42%
Target Name	AChR inhibitor
Calculated MW	368.27
Formulation	C17H21NO3HBr
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CN1CC[C@@]23C=C[C@@H](C[C@@H]2OC4=C(C=CC(=C34)C1)OC)O.BrPathways:NeuroscienceReceptor:AChEBoling

pt:Melting pt:256CSolubility:Water: 36 mg/mL (97.75 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Thomsen T, et al. Eur J Clin Chem Clin Biochem, 1991, 29(8), 487-492.

Note: This product is for in vitro research use only