

Sulfacetamide sodium salt hydrate

Catalog No: #S0087


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Sulfacetamide sodium salt hydrate |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Autophagy inhibitor |
| Calculated MW | 254.2 |
| Formulation | C ₈ H ₉ N ₂ NaO ₃ SH ₂ O |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Immunology
 SMILES: O.[Na].CC(=O)NS(=O)(=O)C1=CC=C(N)C=C1
 Pathways: Autophagy Receptor: DHPS
 Boiling pt: Melting pt: Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hevener KE, et al. J Chem Inf Model. 2009 Jan 27.

Note: This product is for in vitro research use only