

## Pazopanib

Catalog No: #S0097L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pazopanib
Brief Description	Inhibitors
Purification	98.78%
Target Name	c-Kit inhibitor; PDGFR inhibitor; VEGFR inhibitor
Calculated MW	437.52
Formulation	C21H23N7O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CN(C1=CC2=NN(C)C(C)=C2C=C1)C1=NC(NC2=CC(=C(C)C=C2)S(N)(=O)=O)=NC=C1  
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: c-Kit; PDGFR; VEGFR1; VEGFR2; VEGFR3  
 Boiling pt: Melting pt: Solubility: DMSO: 87 mg/mL warmed (198.84 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Harris PA, et al. J Med Chem. 2008, 51(15), 4632-4640.

Note: This product is for in vitro research use only