

Ipratropium bromide monohydrate

Catalog No: #S0098


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ipratropium bromide monohydrate
Brief Description	Inhibitors
Purification	99.37%
Target Name	AChR antagonist
Calculated MW	430.38
Formulation	C ₂₀ H ₃₂ BrNO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory

systemSMILES:CC(C)[N+]1(C)C2CCC1CC(C2)C(=O)C(CO)C1=CC=CC=C1Pathways:NeuroscienceReceptor:mAChRBoling pt:Melting

pt:Solubility:DMSO, Ethanol and Water: 200 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wellington K. Treat Respir Med. 2005; 4(3):215-20; discussion 221-2.

Note: This product is for in vitro research use only