

Atazanavir

Catalog No: #S0100L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Atazanavir
Brief Description	Inhibitors
Purification	98.00%
Target Name	HIV Protease inhibitor
Calculated MW	704.87
Formulation	C38H52N6O7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(C)(C)[C@@H](C(=O)N[C@@H](CC1=CC=CC=C1)[C@H](CN(CC2=CC=C(C=C2)C3=CC=CC=N3)NC(=O)[C@H](C(C)C)NC(=O)OC)O)NC(=O)OC Pathways: Microbiology&virology; Proteases/Proteasome Receptor: HIV Protease Boiling pt: Melting pt: 207-209°C Solubility: DMSO: > 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Robinson BS, et al. Antimicrob Agents Chem Other, 2000, 44(8), 2093-2099.

Note: This product is for in vitro research use only