

Halcinonide

Catalog No: #S0107


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Halcinonide
Brief Description	Inhibitors
Purification	99.00%
Target Name	GR
Calculated MW	454.96
Formulation	C24H32ClFO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC(=O)C=C1CC[C@@H]3[C@@]2([C@H](C[C@]4([C@H]3C[C@@H]5[C@]4(OC(O5)(C)C)(=O)CCl)C)O)FPathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:Solubility:Ethanol: 4 mg/mL warmed (8.79 mM); DMSO: 90 mg/mL warmed (197.81 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhang YD, et al. Am J Surg, 2002, 184(1), 74-77.

Note: This product is for in vitro research use only