

Molsidomine

Catalog No: #S0109


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Molsidomine
Brief Description	Inhibitors
Purification	98.00%
Target Name	NOD
Calculated MW	242.23
Formulation	C9H14N4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular SystemSMILES:CCO/C(=N/C1=C[N+](=NO1)N2CCOCC2)/[O-]Pathways:NF-kbReceptor:NODBoiling pt:Melting pt:140-141CSolubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rosenkranz B, et al. Clin Pharmacokinet. 1996 May;30(5):372-84.

Note: This product is for in vitro research use only