

Idazoxan hydrochloride

Catalog No: #S0115


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Idazoxan hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor antagonist
Calculated MW	240.69
Formulation	C ₁₁ H ₁₂ N ₂ O ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: Cl.C1CN=C(N1)C1COC2=CC=CC=C2O1
 Pathways: Neuroscience; GPCR/G Protein Receptor: Adrenergic Receptor
 Boiling pt: 397°C
 Melting pt: Solubility: Water: 100 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bousquet P, et al. Ann N Y Acad Sci. 1999 Jun 21;881:272-8.

Note: This product is for in vitro research use only