

## Suplatast tosilate

Catalog No: #S0117


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Suplatast tosilate
Brief Description	Inhibitors
Purification	99.00%
Target Name	IL Receptor inhibitor
Calculated MW	499.64
Formulation	C23H33NO7S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:CCOCC(O)COC1=CC=C(NC(=O)CC[S+](C)C)C=C1Pathways:Immunology/InflammationReceptor:Th2 cytokineBoiling pt:Melting pt:84-87CSolubility:Ethanol: 100 mg/mL (200.14 mM); Water: 100 mg/mL (200.14 mM); DMSO: 100 mg/mL (200.14 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mimura T, et al. J Pharmacol Exp Ther, 2005, 314(1), 293-301.

Note: This product is for in vitro research use only