

## Bengenin

Catalog No: #S0122


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bengenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	IL Receptor inhibitor; NF-κB inhibitor; NOD agonist; TNF inhibitor
Calculated MW	328.27
Formulation	C14H16O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory

systemSMILES:COC1=C(C=C2C(=C1O)[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)OC2=O)O)Pathways:Apoptosis;

Immunology/Inflammation; NF-κB Receptor:IL-1β; IL-6; NF-κB; NOD; TNF-α  
 Boling pt:658.9C Melting pt:237-240C Solubility:DMSO: 66 mg/mL (201.05 mM)  
 Appearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Gao XJ, et al. Inflammation. 2015 Jun;38(3):1142-50.

Note: This product is for in vitro research use only