

Selamectin

Catalog No: #S0124


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Selamectin
Brief Description	Inhibitors
Purification	98.00%
Target Name	GluR antagonist
Calculated MW	769.96
Formulation	C43H63NO11
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C[C@H]1CC[C@]2(C[C@@H]3C[C@H](O2)C=C/[C@H]([C@H]/C=C/C=C/4\CO[C@H]5[C@@]4([C@@H](C=C/C5=N/O)C)(=O)O3)O)C[C@H]6C[C@@H]([C@H]([C@@H](O6)C)O)C)O[C@@H]1C7CCCC7 Pathways: Neuroscience Receptor: Glutamate-gated chloride channels Boiling pt: 917°C Melting pt: Solubility: Soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Dryden MW, et al. Parasit Vectors. 2013 Mar 25;6:80.

Note: This product is for in vitro research use only