

Vigabatrin


Catalog No: #S0128

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vigabatrin
Brief Description	Inhibitors
Purification	98.00%
Target Name	ABC antagonist; GABA Receptor agonist
Calculated MW	129.16
Formulation	C6H11NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: NC(CCC(O)=O)C=C
Pathways: Membrane transporter/Ion channel; Neuroscience
Receptor: ABAT; GABA
Boiling pt: 73°C
Melting pt: 209°C
Solubility: Water: 55.1 mg/mL
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Tyacke RJ, et al. *Adv Pharmacol.* 2010;58:373-96.

Note: This product is for in vitro research use only