

Cepharanthine

Catalog No: #S0131


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cepharanthine
Brief Description	Inhibitors
Purification	99.73%
Target Name	Others inhibitor
Calculated MW	606.73
Formulation	C37H38N2O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CN1CCC2=CC3=C(C4=C2[C@@H]1CC5=CC=C(C=C5)OC6=C(C=CC(=C6)C[C@@H]7C8=CC(=C(C=C8CCN7C)OC)O4)OC

OCO3Pathways:OthersReceptor:OthersBoiling pt:Melting pt:140-145CSolubility:DMSO: 100 mg/mL (164.82 mM)Appearance:Yellow

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Huang H, et al. Inflammation. 2014 Feb;37(1):235-46.

Note: This product is for in vitro research use only