

Claudin-3 Polyclonal Antibody

Catalog No: #40750

Package Size: #40750-1 50ul #40750-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Claudin-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Claudin-3 Polyclonal Antibody detects endogenous levels of Claudin-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Claudin-3 around the non-phosphorylation site of Y219.
Target Name	Claudin-3
Other Names	CLDN3; C7orf1; CPETR2; Claudin-3; Clostridium perfringens enterotoxin receptor 2; CPE-R 2; CPE-receptor 2; Rat ventral prostate.1 protein homolog; hRVP1
Accession No.	Swiss-Prot: O15551NCBI Gene ID: 1365
Uniprot	O15551
GeneID	1365;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of COLO cells using Claudin-3 Polyclonal Antibody



Note: This product is for in vitro research use only