

Estradiol valerate

Catalog No: #S0135


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Estradiol valerate
Brief Description	Inhibitors
Purification	99.11%
Target Name	Androgen Receptor agonist
Calculated MW	356.51
Formulation	C23H32O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CCCCC(=O)O[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CCCC4=C3C=CC(=C4)O)CPathways:Endocrinology/HormonesReceptor:Androgen ReceptorBoiling pt:486.2C?Melting pt:144CSolubility:DMSO: 71 mg/mL (199.15 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Takabe K, et al. J Biol Chem, 2010, 285(14), 10477-10486.

Note: This product is for in vitro research use only