

Cefpiramide acid

Catalog No: #S0138


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefpiramide acid
Brief Description	Inhibitors
Purification	98.83%
Target Name	Antibacterial antagonist
Calculated MW	612.64
Formulation	C25H24N8O7S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=CC(=O)C(=CN1)C(=O)N[C@H](C2=CC=C(C=C2)O)C(=O)N[C@H]3[C@@H]4N(C3=O)C(=C(CS4)CSC5=NN=NN5)C(=O)O Pathways: Microbiology & virology Receptor: PBPs Boiling pt: Melting pt: 213-215°C Solubility: Water: insoluble; DMSO: insoluble Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fass RJ, Helsel VL. Diagn Microbiol Infect Dis. 1986 Jan;4(1):77-81.

Note: This product is for in vitro research use only