

Toremifene

Catalog No: #S0139


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toremifene
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor modulator
Calculated MW	405.96
Formulation	C26H28ClNO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: CN(C)CCOC1=CC=C(C=C1)/C(=C(/CCCl)\C2=CC=CC=C2)/C3=CC=CC=C3 Pathways: Endocrinology/Hormones Receptor: ERBo
 ling pt: Melting pt: 108-110°C Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Locci P, et al. J Interferon Cytokine Res. 2001 Nov;21(11):961-70.

Note: This product is for in vitro research use only