

Brinzolamide

Catalog No: #S0142


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Brinzolamide
Brief Description	Inhibitors
Purification	99.50%
Target Name	Carbonic Anhydrase inhibitor
Calculated MW	383.51
Formulation	C ₁₂ H ₂₁ N ₃ O ₅ S ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: CCN[C@H]1CN(S(=O)(=O)C2=C1C=C(S2)S(=O)(=O)N)CCCO Pathways: Metabolism Receptor: Carbonic anhydrase
 IIBoling pt: 586C? Melting pt: Solubility: DMSO: 77 mg/mL (200.77 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Stams T, et al. Protein Sci. 1998 Mar;7(3):556-63.

Note: This product is for in vitro research use only