

Ftaxilide

Catalog No: #S0144


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ftaxilide
Brief Description	Inhibitors
Purification	99.78%
Target Name	Others
Calculated MW	269.3
Formulation	C16H15NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CC1=CC=CC(C)=C1NC(=O)C1=CC=CC=C1C(O)=O
Pathways: Others
Receptor: Others
Boiling pt: 375.3°C
Melting pt: 100°C
Solubility: DMSO: 40mg/mL
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only