

## Calcium folinate

Catalog No: #S0148


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Calcium folinate
Brief Description	Inhibitors
Purification	99.50%
Target Name	DNA/RNA Synthesis inhibitor; Vitamin agonist
Calculated MW	511.51
Formulation	C <sub>20</sub> H <sub>21</sub> CaN <sub>7</sub> O <sub>7</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

SystemSMILES: C1C(N(C2=C(N1)NC(=NC2=O)N)C=O)CNC3=CC=C(C=C3)C(=O)N[C@@H](CCC(=O)[O-])C(=O)[O-].[Ca+2]Pathways: DNA

Damage/DNA Repair; MetabolismReceptor: Thymidylate synthase; VBBoling pt: Melting pt: 240-250CSolubility: Soluble in WaterAppearance: White

PowderRemark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Aschele C, et al. Ann Oncol. 2002 Dec; 13(12):1882-92.

Note: This product is for in vitro research use only