

## Folinic Acid Calcium Salt Pentahydrate

Catalog No: #S0148L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Folinic Acid Calcium Salt Pentahydrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	601.58
Formulation	C <sub>20</sub> H <sub>21</sub> CaN <sub>7</sub> O <sub>7</sub> ·5H <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: NC1=NC(=O)C2=C(NCC(CNC3=CC=C(C=C3)C(=O)NC(CCC(O)=O)C(O)=O)N2C=O)N1 Pathways: Others Receptor: Folate

Analogue Boiling pt: Melting pt: 240-250°C Solubility: Water: 40 mg/mL warmed (66.49 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Van der Wilt CL, et al. Y Res, 1992, 52(18), 4922-4928.

Note: This product is for in vitro research use only