

Eflorate

Catalog No: #S0151


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eflorate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor antagonist
Calculated MW	324.33
Formulation	C19H16O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular SystemSMILES:CCOC(=O)COC1=CC2=C(C=C1)C(=O)C=C(O2)C3=CC=CC=C3Pathways:GPCR/G

ProteinReceptor:β-adrenergic receptorBoiling pt:486.5CMelting pt:Solubility:10 mM in DMSOAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Veronese M, Barzaghi D. Boll Chim Farm. 1985 Mar; 124(3):124-30.

Note: This product is for in vitro research use only