

## Urapidil

Catalog No: #S0165


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Urapidil
Brief Description	Inhibitors
Purification	99.50%
Target Name	5-HT Receptor antagonist; Adrenergic Receptor antagonist
Calculated MW	387.48
Formulation	C <sub>20</sub> H <sub>29</sub> N <sub>5</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: ClCOC1=CC=CC=C1N1CCN(CCCNC2=CC(=O)N(C)C(=O)N2C)CC1  
 Pathways: GPCR/G Protein;  
 Neuroscience Receptor: 5-HT;  $\alpha$ -adrenergic receptor  
 Boiling pt: 549°C  
 Melting pt: Solubility: Soluble in DMSO  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ramage AG. Br J Pharmacol. 1991 Apr;102(4):998-1002.

Note: This product is for in vitro research use only