

Pancuronium bromide

Catalog No: #S0166


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pancuronium bromide
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	732.67
Formulation	C35H60Br2N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC(=O)O[C@H]1C[C@@H]2CC[C@@H]3[C@@H]([C@]2(C[C@@H]1[N+]4(CCCCC4)C)C)C[C@]5([C@H]3C[C@@H]([C@@H]5OC(=O)C)[N+]6(CCCCC6)C).[Br-].[Br-] Pathways: Neuroscience Receptor: AChR; mAChR Boling pt: Melting pt: 214 - 217°C Solubility: Very soluble in water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Stovner J, et al. Br J Anaesth. 1975 Sep; 47(9):949-54.

Note: This product is for in vitro research use only