

Depofemin

Catalog No: #S0168


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Depofemin
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	396.58
Formulation	C26H36O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@@H]2OC(=O)CCC4CCCC4)CCC5=C3C=CC(=C5)OPathways:Endocrinology/HormonesReceptor:ERBoiling pt:532.8CMelting pt:300CSolubility:Ethanol: 55 mg/mL (138.69 mM); DMSO: 79 mg/mL (199.21 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Oriowo MA, et al. Contraception, 1980, 21(4), 415-424.

Note: This product is for in vitro research use only