

Fosfomycin calcium

Catalog No: #S0170


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fosfomycin calcium
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	176.12
Formulation	C ₃ H ₅ CaO ₄ P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: C[C@H]1[C@H](O1)P(=O)([O-])[O-].[Ca+2]
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Water: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Matsumoto T, et al. J Infect ChemOthers. 2011 Feb;17(1):80-6.

Note: This product is for in vitro research use only